



Best™ Integrated Brachytherapy Solutions




Best™ HDR Remote Afterloader

BestMick Applicator



BestMick 10/15/20/25 Cartridges



Best™ Kobold Tandem & Ring Applicator



Best™ Kobold Fletcher Tandem & Ovoid Applicator



Best™ Kobold Henschke Tandem & Ovoid Applicator



Best™ Compact Sonalis™ Ultrasound Imaging System



Best™ Treatment Planning System



Best™ Stepper/Stabilizer

Best™ Esophageal Brachytherapy Applicator



Best™ Breast Double-Balloon Brachytherapy Applicator

Best™ Flexi Needles



15G Flexi Needle



16G Flexi Needle




17G Flexi Needle



15G Flexi Needle with SS Sleeve (Click Fit)



Flexi Needle with Female Hub



Best™ X-ray Marker Cables

Best™ vascular




Vascular Brachytherapy for In-stent Restenosis

VBT Designed for the Cath Lab



- Convenient & Minimal Dosimetry Calculations
- Portable & Reusable
- Compatible with 6F Guide & 7F Guide-Extending Catheters
- Short Treatment Time



Available in 30 mm, 40 mm & 60 mm Source Train Lengths


Best™ Integrated Brachytherapy Solutions

Why Medical Physicists, Radiation Oncologists and Urologists Prefer Best™

Go for the Best™ for your patients!

“A unique aspect of the Best™ seed is its 5 mm length. Best also provides 5 mm laser-cut spacers, which do not have flared edges. With the Best® seed/spacer combination, we can be more confident of accurate seed placement.”¹

Best™ Fiducial Markers



“The source strengths of Best™ seeds are uniformly distributed and the values always fall within the $\pm 3\%$ criterion.”³

“Best™ Model 2301 characteristic x-rays have no significant effect on the photon spectrum and the source emits an essentially pure I-125 decay scheme spectrum.”^{2†}

Best™ Localization Needles with I-125 Seeds



“Best™ Model 2301 seeds are of superior design as they show up more distinctly in transverse ultrasound and fluoroscopic images.”³

“Best™ Model 2301 is the most isotropic of all the I-125 sources currently available. The source is a double encapsulated source.”²

Best™ Palladium-103 Seed



Best™ Iodine-125 Seed



“Delivery of the specific Sk for any given date, delivered on the date requested prior to the procedure. Extremely high level of flexibility, e.g. seed date changes, activity changes, etc. Prompt, courteous customer service, and low price.”⁴

¹ Communication from J. Limmer, UW Cancer Center Wausau Hospital, Wausau, Wisconsin – Dated December 14, 2000. ² B.H. Heintz, R.E. Wallace & J.M. Hevezi, Med. Phys. 28(4) 671-682 (2001).

³ Communication from Khai Lai, Kreiton Oncology Services Shoreline, Washington – Dated March 8, 2001. ⁴ Jeff Limmer, Medical Physicist.



From 0.05 mCi to 10 mCi per I-125 seed, orders are custom prepared.

We ship within **24 hours**, **7 days a week**, sterile and non-sterile, **anywhere** in the world!

Best Medical is the **only** company that makes custom seeds and strands to your exact specifications/prescription.

Best™ Seeds have the best visualization in ultrasound, fluoro, x-ray, CT and MRI, facilitating real-time dosimetry!

